

# Krishna International

## Product Details

---




### ICU Bed Electric KMMEI-101

**Price:** Rs. 0

### **Attributes**

## Product Description

Advance linear actuator system (6 motors, control box, handset, nurse control ACP), extra long 4 pcs side rails with inclination indicators, heavy ABS head/foot panel, 4 section top, 6" central locking castors, SS telescopic IV rod (4 location), sitting position adjustable, leg section extendable, LED reading lamp. Size: 1970-2120L x 900W x 460-740H mm

  
**Krishna**  
— INTERNATIONAL —  
--A Complete Medical Solution--