

Krishna International

Product Details



ICU Bed Mechanical Pediatric

Price: Rs. 0

Attributes



Product Description

4 section mattress area (CRCA sheet), detachable SS wooden laminated head/foot bows, collapsible/drop-down safety side rails, backrest/knee rest/height adjustment/trendelenburg via smooth crank mechanism with folding SS handles, corner PU buffers, heavy duty IV pole (4 locations), single wheel diagonal locking castors (125mm). Overall: 1600L x 850W mm

Krishna
— INTERNATIONAL —
--A Complete Medical Solution--