

Krishna International

Product Details



ICU Bed Mechanical Five Function (KMMEI-105)

Price: Rs. 0

Attributes

Product Description

4 section mattress area (CRCA sheet), detachable ABS head/foot bows, 4 pcs individual ABS side rails with angle indicator, corner PU buffers, saline stand provision (4 locations), 125mm diagonal locking castors, back rest/knee rest/trendelenburg/reverse trendelenburg maneuverable via 4 separate crank mechanisms from foot-end. Frame size: 2030L x 900W x 600-800H mm

