

Krishna International

Product Details




ICU Bed Electric KMMEI-101

Price: Rs. 0

Attributes

Product Description

Advance linear actuator system (6 motors, control box, handset, nurse control ACP), extra long 4 pcs side rails with inclination indicators, heavy ABS head/foot panel, 4 section top, 6" central locking castors, SS telescopic IV rod (4 location), sitting position adjustable, leg section extendable, LED reading lamp. Size: 1970-2120L x 900W x 460-740H mm


Krishna
— INTERNATIONAL —
--A Complete Medical Solution--